

Circulaire® Touchscreen Fume Cupboards



CT 1100

CT PRO 1400

The Circulaire® Touchscreen Filtration Fume Cupboards solutions feature the very latest in air filtration technology.

Fully compliant with COSHH regulations and international standards including BS7989:2001 for Filtration Fume/Particulate Cupboards.

Monmouth Circulaire® Fume Cupboards use as low as 100 watts of power to ensure low energy costs. The units employ extra-large capacity Activated Carbon filters for removal of fumes and/or HEPA filters for particulate removal. The filters are up to 30% bigger than those of rival Fume Cupboards and are significantly deeper making fume removal more efficient and effective.

Our revolutionary Visionaire® Touchscreen Control Interface is integrated into the cabinet to control & monitor all aspects of product operation.



CARBON FILTER

Size Specification (W x D x H)

CT 800	CT 1100	CT 1400	CT 1800
800mm x 700mm x 1284mm	1100mm x 700mm x 1284mm	1400mm x 700mm x 1284mm	1800mm x 700mm x 1284mm
CT PRO 800	CT PRO 1100	CT PRO 1400	CT PRO 1800
800mm x 700mm x 1345mm	1100mm x 700mm x 1345mm	1400mm x 700mm x 1345mm	1800mm x 700mm x 1345mm