

Circulaire® Fume Cupboards



C 650

**BRITISH
MANUFACTURED**

**FIVE YEAR
WARRANTY**

**ULTRA DEEP
CARBON
FILTERS**

**CUSTOMISABLE
FOR YOUR
APPLICATION**

The Circulaire® Filtration Fume Cupboard is fitted with ultra-deep Activated Carbon or HEPA filters to guarantee maximum fume filtration efficiency.

All Monmouth products are designed to have both low noise levels and low energy consumption to help your company stay carbon neutral. LED Lighting further improve efficiency.

Fully compliant with COSHH regulations and international standards including BS7989:2001 for Filtration Fume/Particulate Cupboards.

To ensure the cabinet is functioning correctly, every model has an animated colour LCD display showing the airflow is working within safe and calibrated parameters.

Size Specification (W x D x H)

C 550	C 650	C 900
550mm x 600mm x 1121mm	650mm x 600mm x 1121mm	900mm x 600mm x 1121mm