


# Ductaire® Ducted Fume Cupboards



**BRITISH  
MANUFACTURED**

**FIVE YEAR  
WARRANTY**

**FULLY  
CHEMICALLY  
RESISTANT**

**CUSTOMISABLE  
FOR YOUR  
APPLICATION**

**The Ductaire® Series has been designed with the end-user in mind and can be installed into virtually any laboratory.**

Compliant with BS EN 14175, the Ductaire® offers excellent value for money and a variety of features to provide an excellent Ducted Fume Cupboard.

The Fume Cupboard has a high chemical resistance and a fully welded polypropylene construction, offering an integral spillage containment area. Each unit is fully tested and certified and offers exceptional airflow characteristics and user safety. The front glazing folds up to allow full access to the working chamber and all electronics are fully isolated from this chamber.

## Size Specification (W x D x H)

D 700	D1000	D 1200	D 1500	D 2000	D 3000
700mm x 700mm x 970mm	1000mm x 700mm x 970mm	1200mm x 700mm x 970mm	1500mm x 700mm x 970mm	2000mm x 700mm x 970mm	3000mm x 700mm x 970mm