

# Ductaire® Titan Ducted Fume Cupboards



DT1000

**The Ductaire® Titan offer the best possible features and quality in a Ducted Fume Cupboard.**

Applications include; Laboratory, Pharmaceutical and Industrial Analytical Research as well as Industrial, Medical and Educational Establishments requiring extraction systems for operator protection.

## Standard Features Include:

- Outer Carcass fabricated from Zintec Mild Steel, Powder Coated in Textured Matt White (RAL9016).
- 6mm Toughened Glass Sash, includes stops set at 500mm and 25mm.
- Alarm System - LCD Air Speed read-out, Fan stop/start and programmable relay outputs and inputs.
- Interior 55W, Low Energy Lighting.
- Removable back baffle for easy cleaning and maintenance.
- Fully BS EN 14175 & COSHH compliant.

Their robust design and quality fabrication mean the systems offer years of uninterrupted use and like all Monmouth products, the range comes with a 5-Year Warranty.

## Size Specification (W x D x H)

DT 1000	DT 1200	DT 1500	DT 1800	DT 2000
1000mm x 900mm x 2400mm	1200mm x 900mm x 2400mm	1500mm x 900mm x 2400mm	1800mm x 900mm x 2400mm	2000mm x 900mm x 2400mm